



**HAL**  
open science

## Identification of plastic additives: Py/TD-GC-HRMS method development and application on food containers

Fleurine Akoueson, Chaza Chbib, Armance Brémard, Sébastien Monchy, Ika Paul-Pont, Périne Doyen, Alexandre Dehaut, Guillaume Duflos

### ► To cite this version:

Fleurine Akoueson, Chaza Chbib, Armance Brémard, Sébastien Monchy, Ika Paul-Pont, et al.. Identification of plastic additives: Py/TD-GC-HRMS method development and application on food containers. *Journal of Analytical and Applied Pyrolysis*, 2022, 168, pp.105745. 10.1016/j.jaap.2022.105745 . anses-03846369

**HAL Id: anses-03846369**

**<https://anses.hal.science/anses-03846369>**

Submitted on 10 Nov 2022

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Copyright

1 **Identification of plastic additives: Py/TD-GC-HRMS method**  
2 **development and application on food containers.**

3 Fleurine Akoueson<sup>1,2</sup>, Chaza Chbib<sup>3</sup>, Armance Brémard<sup>1,5</sup>, Sébastien Monchy<sup>3</sup>, Ika  
4 Paul-Pont<sup>4</sup>, Périne Doyen<sup>2</sup>, Alexandre Dehaut<sup>1</sup>, Guillaume Duflos<sup>1</sup>.

5 <sup>1</sup> ANSES – LSAI, Boulevard du Bassin Napoléon, 62200 Boulogne-sur-Mer, France.

6 <sup>2</sup> Univ. Littoral Côte d'Opale, UMR 1158 BioEcoAgro, EA 7394, Institut Charles Viollette, USC  
7 ANSES, INRAe, Univ. Lille, Univ. Artois, Univ. Picardie Jules Verne, Uni. Liège, F-62200,  
8 Boulogne-sur-Mer, France

9 <sup>3</sup> Univ. Littoral Côte d'Opale, CNRS, Univ. Lille, UMR 8187, LOG, F 62930 Wimereux, France

10 <sup>4</sup> Univ Brest, CNRS, IRD, Ifremer, LEMAR, F-29280 Plouzané, France

11 <sup>5</sup> Université Paris Saclay, Bât. 349, 91405 Orsay Cedex, France

12 Corresponding author: [Guillaume.duflos@anses.fr](mailto:Guillaume.duflos@anses.fr)

## Highlights

- A single method for the analysis of OPAs using Py/TD-GC-HRMS has been implemented
- The developed method allowed to create a plastic additive database
- 49 additives molecules out of the selected 57 have been successfully identified
- 10 plastic additives were found in PP and PLA food packaging

13 **Abstract**

14 The interest of studying organic plastic additives (OPAs) is growing since many of them  
15 are considered as toxic for human beings, marine organisms and have been found in  
16 the ecosystem. The impact and the extent of this pollution is not fully understood  
17 because of the lack of methods able to comprehensively analyze the wide diversity of  
18 OPAs. This study aims at offering a quick and efficient way to identify specific harmful  
19 plastic additives in polymer materials. For this purpose, a general method, for  
20 identifying 49 OPAs, has been successfully developed using pyrolyzer/thermal  
21 desorption-gas chromatography coupled to a high-resolution mass spectrometer  
22 (Py/TD-GC-HRMS). The retention indices were calculated for each compound and a  
23 high-resolution mass spectra database for plastic additives was implemented. This  
24 method was then applied to characterize the OPAs included in food containers that are  
25 frequently found in the environment. Four samples made of polypropylene (PP) and  
26 polylactic acid (PLA) coming from two distinct suppliers were then selected. A total of  
27 10 additives, which were plasticizers (phthalates and adipates), antioxidants  
28 (bisphenol B) and phosphorous flames retardants, were found in these food packaging.  
29 Generally, higher peak signals were detected in the PP samples in comparison to PLA  
30 samples, except for tributyl phosphate and bisphenol B that were exclusively found in  
31 PLA samples. Moreover, the experimentations pointed out additives differences  
32 between suppliers, for products intended for the same use and made of the same  
33 polymer, for both PP and PLA samples highlighting the complexity of identifying and  
34 generalizing the additive content of plastic materials.

35 **Keywords:** Phosphorous flame-retardants, Plasticizers, antioxidants, UV stabilizers,  
36 Py/TD-GC-HRMS, analytical method development.

## 1. Introduction

Plastic pollution is a threat for the environment and marine organisms [1]. The presence of plastics is ubiquitous in the environment. In 2020, plastic production reached 307 million metric tons worldwide [2], of which 24% were used for food and beverages packaging [3]. These items are in the top ten of the most common plastic waste found in the environment [4].

A key concern of plastic pollution is linked to the risks to both ecosystem and human health. These adverse effects classified in several categories including the “physical” effects related to the inhalation and ingestion of plastic debris and the “chemical” effects related to the harmful chemicals associated with plastics such as the organic plastic additives (OPAs) included in the polymer [5,6]. While the “physical” impacts of plastics are well described in the literature, the “chemical” impacts, partly related to the presence of OPAs, is less documented and limited information is available for the majority of them. Several categories and families of OPAs exist, sorted by decreasing total amounts in polymers: plasticizers, flame retardants, colorants, antioxidants, heat stabilizers, UV stabilizers or bio-stabilizers (biocides), lubricants, foaming agents, and antistatic agents. These chemicals are added to polymers during their formulation, without being chemically bound to polymers [7], in order to confer them specific properties and advantages. However, there are increasing concerns about their impacts on marine organisms or human health since some of them were found to be toxic [7]. Indeed some phthalates, nonylphenols (NPs), Bisphenols, as well as some brominated flames retardants (BFRs) have been identified as carcinogens, mutagens and may have detrimental impacts on the reproduction being endocrine disrupting chemicals [1,8-11]. As most of the OPAs are not chemically bound to the polymeric matrix, they can leach into the surrounding environment or can be transferred to marine organisms or humans upon plastic debris uptake [5,12-14]. Additionally, OPAs may migrate from the packaging leading to contamination of the food packed [15-18]. As a result, the use of OPAs in plastics has become controversial and the assessment of their nature and toxicity is of utmost importance.

A considerable number (more than 5300) of polymer formulations are available on the market, and more than 4000 known chemicals are used in the formulation of plastic packaging, i.e. intentionally added substances (IAS). Additionally, plastic packaging may contain non-intentionally added substances (NIAS) [19,20]. That demonstrates the chemical complexity of plastic polymers. The identity of most IAS and NIAS are unknown making difficult to assess the safety of plastic materials, especially in regards with food contact materials (FCM). In order to broaden understanding of the chemical impact of these additives, it is necessary to have a more-detailed assessment of plastic packaging associated chemicals. This involves the development of robust analytical method to identify and quantify them. The analysis of OPAs in polymers is an analytical challenge because of the wide array of their chemical compositions

73 and the diversity of the molecules mixed within a polymer. Also, the most frequently used  
74 processes, e.g. solvent extraction (SE), present constraints. These methods are time  
75 consuming for the sample preparation, generate wastes and involve the use of toxic solvents.  
76 Besides, there is no universal solvent able to dissolve all types of polymers [21-24]. Under  
77 such circumstances, the need of avoiding the complexity of analytical processes is of concern.  
78 Pyrolysis coupled to a gas chromatography and a mass spectrometry (Py-GC-MS) is an  
79 analytical tool that enables relatively quick and easy analyses [25]. Moreover, this technique  
80 minimizes the sample preparation steps, limiting the potential contamination prior to the  
81 analysis and does not require the use of toxic solvents [21,26-28]. Analytical pyrolysis also  
82 offers a major advantage. It permit to successively analyze both the OPAs and the polymer  
83 with a multi-step approach including an initial thermal desorption (TD) step at sub pyrolytic  
84 temperatures (i.e. approximately <500 °C) to characterize the additives regardless of the  
85 polymer, and a second pyrolysis step at higher temperatures (i.e. >500 °C) to analyze the  
86 polymeric matrix [29-31]. Py-GC-MS has already been reported as beneficial to thermally  
87 desorb and screen plastic additives in several studies using pyrolyzer/thermal desorption  
88 (Py/TD) GC-MS methods [32-34] and is a recognized method for the fast identification of  
89 OPAs. Akoueson et al., (2021) review, listed a wide array of OPAs studied using Py and Py/TD-  
90 GC-MS [25]. However all the developed methods focused on few analytes from a type of  
91 additives. No general method exists to analyze a wider diversity of additives. Thus, we aimed  
92 to develop a TD method able to identify diverse harmful, emergent, or widely used OPAs.

93 The objectives of the present study are: (i) set up a high-resolution Mass spectral library which  
94 includes a selection of OPAs used as plasticizers, flame retardants, antioxidants and  
95 stabilizers, (ii) implement a general analytical method using Py/TD-GC-HRMS to identify and  
96 detect the selected plastic additives, and (iii) apply the method to several items of food  
97 packaging materials in order to characterize their OPAs content.

## 98 **2. Materials and methods**

### 99 **2.1. Implementation of a general GC-HRMS method for OPAs** 100 **detection**

#### 101 **2.1.1. Targeted additives, reagents and preparation of the samples**

102 According to different criteria, detailed in part 3.1, 56 additives were selected. The additives  
103 and their corresponding suppliers are listed in Table 1. A standard mixed solution of C<sub>3</sub>–C<sub>33</sub> n-  
104 alkanes (Cat. No. 40147-U) was purchased from RESTEK (Lisses, France). Analytical GC  
105 grade solvents were used such as methanol, acetone, toluene or acetonitrile and were all  
106 purchased from Sigma-Aldrich Co. (Saint-Quentin-Fallavier, France).

107 Diluted solutions of each additives stock solutions were prepared in suitable solvents at  
108 concentrations ranging from 0.05 to 100 mg.L<sup>-1</sup>. The diluted solutions were then transferred  
109 and conserved in 2 mL GC autosampler vials with a PTFE coated seal. An aliquot of 1 µL of  
110 each additive solution was injected in the GC-HRMS using a microliter syringe at increasing  
111 concentrations (0.12, 0.25, 0.5, 1, 2.5 and 5 mg.L<sup>-1</sup>) until visualization of the limit of detection  
112 additives when the signal to noise (s/n) was > 3. For the mixed n-alkanes, 1 µL of the 10 mg.L<sup>-1</sup>  
113 solution was also injected.

### 114 **2.1.2. GC-HRMS method and parameter settings**

115 Samples were analysed using a GC Trace 1310-MS Orbitrap Q-exactive from ThermoFisher  
116 Scientific (Les Ulis, France) equipped with a TriPlus RSH auto-sampler. The compounds in  
117 standard solutions were separated on a 30 m × 0.25 mm × 0.25 µm RXi-5ms capillary column  
118 acquired from Restek (Lisses, France). The injection volume was 1 µL.

119 Based on the literature, several methods (Supplementary data S1) were compiled in order to  
120 propose a compromise, i.e. a general method suitable for the 56 selected additives, being  
121 aware of its limitations regarding the analysis sensitivity for each additive.

122 Analytical conditions are detailed in Table 2.

123 Chromatograms were analyzed using Thermo Scientific Xcalibur software (4.3.73.11) and with  
124 the help of an automated search for the markers by deconvolution using TraceFinder™  
125 software (5.0.889.0). The identity of the additives was confirmed using target and confirming  
126 ions obtained from the literature [35-40] and the National Institute of Standards and Technology  
127 (NIST) Database.

### 128 **2.1.3. Construction of the database**

129 Following the acquisition, and after obtaining robust chromatograms signals for all the  
130 analytes, a database containing mass spectra for all the detected additives was implemented  
131 using TraceFinder software including quantitative and qualitative ions for each molecule, as  
132 specific as possible. The HRMS files are available on Zenodo [41].

133 Following the analysis of the standard additives solutions and the mixed *n*-alkane solution, the  
134 retention time (RT) of the additive analytes and the *n*-alkanes were used to calculate the  
135 Kováts retention indices (RI) following equation (1). Respective RI were calculated [42] and  
136 registered in the database. The RI were validated when the calculated values were in the range  
137 of RI values found in the literature.

$$138 \quad RI = 100 \times (C_{n+1} - C_n) \times \left[ \frac{(t_x - t_n)}{(t_{n+1} - t_n)} \right] + 100 \times C_n \quad (1)$$

139 With the carbon numbers of *n*-alkanes eluted before ( $C_n$ ) and after ( $C_{n+1}$ ) the corresponding  
140 target molecule, retention time of the target molecule expressed in min ( $t_x$ ) and retention times  
141 for  $C_n$  ( $t_n$ ) and  $C_{n+1}$  ( $t_{n+1}$ ) for the flanking *n*-alkanes expressed in min.

142  
143 The RI and the specific qualitative and quantitative ions are gathered for each additive in  
144 Table 3.

## 145 **2.2. Thermodesorption method development for the detection of** 146 **OPAs using Py/TD-GC-HRMS**

147 A TD step using a Pyrolyser, CDS Analytical Pyroprobe 6150 from Quad Service (Oxford, PA  
148 US), was added to the GC-HRMS developed method.

### 149 **2.2.1. Sample preparation optimization**

150 Four protocols of sample preparation (Supplementary data S2) were previously tested in order  
151 to ease the analytical process and to reduce the high volatility of several low molecular weight  
152 additives. Briefly, the analysis of the selected plastic additives have been done by introducing  
153 1  $\mu\text{L}$  of the additives stock solutions at concentrations ranging from 5 to 100  $\text{mg}\cdot\text{L}^{-1}$ . The  
154 solutions were either added in an empty pyrolysis quartz tube or in a pyrolysis quartz tube  
155 containing three types of different materials used as “solution trapping matrix” such as a piece  
156 of quartz filter, a piece of quartz wool or a small amount of diatomaceous earth (silicium  
157 dioxide). These media were selected for their documented chemically inert features and their  
158 thermal resistance [43,44]. In brief, considering the results (S2), the quartz filter was chosen  
159 as a suitable « trapping matrix » since it permit to stabilize the compounds despite the waiting  
160 time between the sample preparation and the analysis

### 161 **2.2.2. Implementation of a general method using Py/TD-GC-HRMS device**

162 The analyses were carried out using a Pyrolyser CDS Analytical Pyroprobe 6150 from Quad  
163 Service (Oxford, PA US) coupled to a GC-HRMS device (Trace 1310-MS Orbitrap Q-exactive,  
164 ThermoFisher Scientific) thanks to a heated transfer line with the method presented in part  
165 2.1.2.

166 According to the preliminary results obtained following the sample preparation optimization  
167 tests presented part 2.2.1., 1  $\mu\text{L}$  of the additives standard solutions was introduced in a  
168 pyrolysis quartz tube containing a piece of quartz filter. The pyrolysis quartz tube and the quartz  
169 filter were both heated beforehand at 1000 °C during 15 seconds to get rid of any residual  
170 compounds. Sample tubes were then placed in the pyrolyzer auto-sampler and analysis was  
171 performed. Selection of the thermal desorption temperature was carried out on the additives  
172 standard solutions by testing three desorption temperatures: 250 °C, 350 °C and 450 °C, with



173 three replicates. The TD time was fixed at 60 seconds based on the literature, being the most  
174 frequently used TD time (Table S1).

175 To ensure the reliability of the results and the absence of contamination in the analytical device  
176 two precautions were taken for each batch: (i) a 1000 °C auto-clean of the pyrolysis chamber  
177 during 15 seconds between each analysis; (ii) and five empty pyrolysis quartz tube without  
178 sample, hereafter referred as “blank”, were pyrolyzed at the beginning of each run.

## 179 **2.3. Application: OPAs identification in food contact materials**

### 180 **2.3.1. Sample selection and preparation**

181 Two types of food packaging items made of polypropylene (PP) and polylactic acid (PLA),  
182 coming from two different suppliers, tagged with either A or B, have been selected. The  
183 selection of PP was based on its high production and high presence in the marine environment  
184 [3,19]. The PLA was selected for being the most used bio-based alternative in order to compare  
185 its additive content to a petro-based polymer [45].

186 Punches, i.e. cylinders of 1 mm diameter, were cut in the food packaging items using biopsy  
187 punches from Farla-Medical (Antwerpen, Belgium), in order to have a homogeneous sample  
188 taking process. Thickness of the cylinders was measured using a binocular microscope  
189 (Olympus SZX16) and OLYMPYS cellSens Dimention 3.2 software (n=3). The volumes and  
190 the masses ( $\mu\text{g}$ ) of each particle were calculated using the diameters and the thickness  
191 measurements. The cylinders were then introduced in pyrolysis sample cups. Sample tubes  
192 were placed in the pyrolyzer auto-sampler and analysis was performed (n=3).

### 193 **2.3.2. Py/TD-GC-HRMS analyses**

194 In order to thermally desorb the additives contained in the samples, each cylinder was  
195 successively analyzed three times at a thermal desorption temperature of 350 °C to maximize  
196 the extraction of the chemicals present in the samples. An approximation of the extraction  
197 efficiency percentage (%) was calculated dividing the area recorded at the first step by the total  
198 area over the three steps.

199 All analyses were carried out using the method previously developed and presented part 2.1.2.

200 The selected OPAs (Table 1) were identified on the basis of the RT, RI, quantitative and  
201 qualitative ions and ions ratios with the help of the previously implemented database. The  
202 obtained signals were corrected taking in account the procedural blanks, i.e., pyrolysis tubes  
203 without sample, the calculated mass of the samples and the extraction efficiency of the method.

## 2.4. Statistical analysis

All statistical analyses were performed using R-Studio software (1.4.1106) [46]. For all the analysis, the normality and the homoscedasticity of the distribution were verified before carrying out ANOVA follow by Tuckey post Hoc test (adding the car package (3.0-12) [47]). Assuming that one of the hypothesis was not verified, a non-parametric Kruskal-Wallis test was performed. Kruskal Wallis tests were followed by a Nemejni's post hoc test. Post hoc tests were achieved using agricolae (1.3-5) [48] and PMCMR (4.4) [49] packages. Mean differences were considered as significant when p-value < 0.05.

## 3. Results and discussion

### 3.1. Preliminary work and implementation of a HR mass spectra database

A total of 56 additives were selected as relevant compounds. This selection was based on a bibliographic work including previous publications, the European Food Safety Authority's (EFSA) list containing their evaluation on substances for FCM and the European chemical agency's (ECHA) additives list containing 400 substances imported in the European Union that are used in plastic products with a trade volume higher than 100 ton/year (Table 1) [1,7,23,25,26,34,45,50-55]. This bibliographic work was performed based on different priority criteria such as: the use of additives in plastic and plastic packaging, a balance between their toxicity and their concentration in polymers, the environment and marine organisms, and finally, the feasibility of GC-MS analysis on these additives.

A standard solution containing all targeted compounds was prepared to set up a suitable GC–HRMS method for their simultaneous analysis. This preliminary work was carried out using GC-HRMS with the method presented part 2.1.2 and allowed the detection of all the selected additives (Figure 1). The spectra of each additives peaks in figure 1 permitted the implementation of our own plastic additives database [41].

The identification of the additives not already included in NIST database was confirmed by the presence of their respective quantitative and qualitative characteristic ions presented in table 3. An additional confirmation for this set of additives has been made using Kováts RI (Table 3).

In this work, as the analytes were partly selected on their amenability to GC-MS analysis, some additives known to be toxic and widely used in plastic materials (e.g. Tetrabromobisphenol A and the hexabromocyclododecane isomers ( $\alpha$ ,  $\beta$  and  $\gamma$ -HBCD)) were not targeted. Indeed, depending on their chemical nature, certain OPAs may be complex to analyze and are hard to mobilize [34]. These analytical limitations can be overcome with a derivatization step that will

238 modify the sample and product more volatile derivative [56]. It was not employed in the present  
239 study, which consisted to implement a fast and easy method. However, using a derivatization  
240 approach or leading analysis using liquid extraction could be complementary to have a more  
241 complete analysis.

### 242 **3.2. Optimization and characterization of a Py/TD-GC-HRMS method** 243 **for the detection of OPAs**

244 Three temperatures were tested at this step: 250 °C, 350 °C and 450 °C. The pyrolyzer thermal  
245 desorption temperature had a significant impact on the peaks area of several additives. Figure  
246 2 illustrates the relevant temperature patterns for each additive families. The detailed results  
247 for each additive are presented in supplementary data S3.

248 The areas of most flame-retardants (9 out of 14) showed no significant difference between the  
249 three tested temperatures. Areas were significantly higher at 350 °C for three additives (TCEP,  
250 TCPP and TDCPP), whereas the area of the 2,4,6-TBP and BTBPE were significantly higher  
251 respectively at 250 °C and 450 °C. The analysis of the plasticizers additives indicates  
252 significantly higher areas at 350 °C for 13 additives (BBP, DBP, DCHP, DEHP, DEP, DHP,  
253 DIBP, DIDP, DIHP, DINCH, DNOP, DNP, DOA) out of 16, and at 450 °C for the ATBC. No  
254 significant difference has been detected for three additives between the three tested  
255 temperatures (DAP, DIHP and DMP). Concerning the antioxidants, areas were significantly  
256 higher at 450 °C for six additives (NPs, 4-NP, NP10E, 4NP10E and Irganox 1081®), at 350°C  
257 for the Lowinox TBM 8® and at 250 °C for the 4-t-OP. Five additives showed no significant  
258 difference between the tested temperatures. Finally, no significant difference was observed for  
259 the three UV stabilizers.

260 Overall, based on these results, 350 °C seems to be the more relevant thermal desorption  
261 temperature for the analysis of the selected additives, either being significantly optimal or  
262 because no significant differences has been observed between the three temperatures. In line  
263 with our results, other studies generally used a thermal desorption temperature of 350 °C  
264 [25,27,52]. Moreover, the temperature is one of the key factors for a successful extraction and  
265 subsequent detection [27]. Indeed, it affects the degree of fragmentation and the generated  
266 pyrolysis products, its selection must allow a suitable desorption of the additive without altering  
267 it. At too high temperatures, usually above 500 °C, the polymeric matrix begins to degrade and  
268 the excess of polymer fragments interfere with the identification of characteristic additive  
269 fragments, which make the analysis of the results complicated [21] and may lead to  
270 misinterpretation of the results.

271 The Py/TD-GC-HRMS method was used to analyze 56 additives at various concentration  
272 levels (1, 5, 10, 25, 50 and 100 mg.L<sup>-1</sup>). Eight additives, marked with “\*” in the Table 1, were

273 not detected even when the highest concentration was injected in the pyrolyzer. Among them,  
274 three phosphorous flame-retardants (TCP, TCrP and TToP), two plasticizers (bisphenol S and  
275 TEHPA), two UV stabilizers (Uvinul 3049<sup>®</sup> and 3008<sup>®</sup>) and one antioxidant (NP2OE). The 49  
276 additives identified have been detected at 5 ng except for the BDEs, the alkylphenols and the  
277 Irganox 1081<sup>®</sup> at 10 ng; for DIDP, DINP, DIHP and the BTBPE at 25 ng; for Lowinox TBM 6<sup>®</sup>  
278 at 50 ng; and for the Irganox 1076<sup>®</sup> with a detection threshold at 100 ng (Table 1).

279 A previous study [53] mentioned the complicated analysis of phosphorous flame-retardants.  
280 Indeed, some analytes are more difficult to analyze than others using Py/TD-GC-HRMS  
281 according to their physico-chemical properties [25]. For instance, TEHPA and Uvinul 3008<sup>®</sup>  
282 start to decompose at a temperature of 250°C and 350°C respectively, unlike BPS that has a  
283 low volatility and is more thermally stable than its bisphenol's analogs [57].

284 Increasing the thermal desorption temperature might help to volatilize the most stable  
285 additives, such as the BPS. Nonetheless, the temperatures selected in the study did not permit  
286 the desorption of the BPS. Other studies managed to analyze this molecule, for instance  
287 Becerra and Odermatt (2012) conducted their analysis with a temperature of 500 °C and added  
288 derivatizing agents using trimethyl ammonium hydroxide (TMAH), to increase the volatility of  
289 the additive [58]. However, in the present study, it was a thoughtful choice not to derivatize the  
290 samples, since the aim was to implement a quick method with minimal sample preparation  
291 steps. Some polymers, such as polystyrene (PS) or polyvinyl chloride (PVC), start to  
292 decompose at relatively low temperatures [21]. These limitations explain the selection of a  
293 method that does not permit the analysis of all additives.

### 294 **3.3. Application of the method on plastic food packaging**

#### 295 3.3.1. Integration of samples mass variability

296 The masses ( $\mu\text{g}$ ) of each sample cylinders were calculated using their diameter and thickness  
297 as well as the density of each polymer type. The average calculated mass of PP-B, i.e. 203  
298  $\mu\text{g}$ , is significantly different ( $p\text{-value} < 0.05$ ) from the three other samples, i.e. 338  $\mu\text{g}$ , 317  $\mu\text{g}$   
299 and 276  $\mu\text{g}$  respectively for PP-A, PLA-B and PLA-A. Moreover, a high standard deviation (sd)  
300 is noticed especially for PLA-B sample and can be explain by thickness variations according  
301 to the location where the sample was cut out on the plastic container (supplementary data S4).  
302 In that respect, the respective masses of each analyzed food packaging cylinders were taken  
303 into account in order to normalize the additives signal measured in each sample.

#### 304 3.3.2. Efficiency characterization

305 The method was developed on additives standard solutions. However, once additives are  
306 included in a polymeric matrix, some interactions can occur between the additives and the  
307 matrix [21]. Hence, the extraction efficiency of the method was directly calculated on the

308 matrices of interest, i.e. following the analysis of PP and PLA food packaging samples, in order  
309 to take into account the potential interactions.

310 Each sample were submitted to three successive heating at a thermal desorption temperature  
311 of 350 °C. No additive was detected during the third analysis, indicating that the entirety of the  
312 additive included in the samples has been desorbed and analyzed during the first or second  
313 pyrolysis step (supplementary data S5). Consequently, the sum of the areas obtained during  
314 the two first analysis steps is considered as representing 100% of the additives amount  
315 included in each sample. The extraction efficiency percentage of each step was therefore  
316 calculated based on this assumption.

317 The results of the extraction efficiency after a first step of desorption were more than 80% for  
318 most compounds (Table 4, bold printed characters), with a few exceptions. Five additives  
319 required a total of two desorption steps to be fully extracted from the samples, e.g. DMP, DEP  
320 (for PP samples), DEHA (PP-A), TBP (PLA-A and PP-A) and TCPP, with average extraction  
321 efficiency percentages ranging from 50 to 80% or <50% upon the first desorption step (Table  
322 4). Carrying out two desorption steps doubles the time of the analysis. For this reason, the  
323 developed method only includes one thermal desorption step at 350 °C. The average  
324 extraction efficiency percentages of each additive for each sample were taken into account in  
325 order to correct the additives signal recorded. However, care must be taken about the  
326 completeness of the extraction with this TD method. Although the punch is thin and the cross-  
327 section expose the core of the material. It can also be assumed that only the additives on the  
328 first layers of the punch could have been desorbed.

### 329 3.3.3. Identification of the selected OPAs in PP and PLA food packaging samples

330 The additives amounts detected in the food packaging samples are expresses in Area  
331 (Arbitrary unit (A.U)) per µg of plastic. The values were corrected with the extraction efficiency  
332 percentage and the calculated mass of each particle analyzed.

333 Ten OPAs of interest have been detected in the selected food packaging samples (Figure 3),  
334 mostly plasticizers (eight molecules) and flames retardants (two molecules). Two plasticizers,  
335 DEP and DIBP and a phosphorous flame retardant, TCPP, have been detected in all samples  
336 regardless of the supplier or polymer type. DMP and DBP plasticizers have been identified in  
337 all the samples except PLA-A, and DCHP was identified in all the samples except PP-B.  
338 Finally, the plasticizer DINP has only been detected in the two PP samples and Bisphenol B  
339 exclusively in the two PLA samples. TBP, a phosphorous flame retardant, has been identified  
340 both in PP-A and PLA-B only.

341 The area of DCHP is equivalent among all samples regardless of the supplier or polymer type.  
342 Overall, peaks areas for a majority of the OPAs are higher for the PP samples than for the PLA

343 (six out of ten; Figure 3). DBP and DEP peak areas are significantly higher in PP-A samples  
344 than in PLA-A and B samples. The areas of DEHA, DIBP, DINP and TCPP in both PP samples  
345 are significantly higher than in PLA samples. On the other hand, the peak area of the  
346 phosphorous flame retardant TBP is significantly higher in the PLA-A than in both PP samples.

347 Looking at the differences between PLA-A and B (Figure 3), the PLA-B contains more additives  
348 than the PLA-A. Nine OPAs were identified in the PLA-B, including seven plasticizers (BPB,  
349 DBP, DCHP, DEHA, DEP, DIBP and DMP) and two flames retardants (TBP and TCPP), while  
350 only seven OPAs were identified in the PLA-A in which the DEHA and DMP have not been  
351 detected. In addition, a trend appears that the areas of the OPAs peaks found in PLA-B are  
352 mostly higher than those found in the PLA-A, except for the DCHP where the areas are  
353 identical. This trend is significant for the Bisphenol B (15 times higher in PLA-B) and the DEHA  
354 (absent in the PLA-A while detected in the PLA-B). The same pattern is observed when looking  
355 at the PP-A and B samples results (Figure 3). The PP-A contains two additives more than the  
356 PP-B. Nine OPAs have been detected in the PP-A, including seven plasticizers (DBP, DCHP,  
357 DEHA, DEP, DIBP, DINP and DMP) and two flames retardants (TBP and TCPP), while the  
358 PP-B sample only contains seven OPAs (DCHP and TBP missing). The peaks areas of the  
359 OPAs identified in the PP-A samples tends to be higher than in the PP-B. It is significant for  
360 the DBP (three times higher in PP-A), the DEHA (8x), the DIBP (x1.5), DMP (x2) and TCPP  
361 (x4).

362 However, regarding DEP and DIBP (respectively for PLA and PP, and for PLA only), it is  
363 difficult to draw conclusions. Indeed, a large spread in the data can be observed as these  
364 additives show a wide 95% confidence interval (IC 95 %). This dispersion can be explained by  
365 a lack of homogeneity in the distribution of additives within the same sample. This would then  
366 require the multiplication of replicates to be able to take these disparities into account.

367 These results lead to a few remarks:

368 First, a petro-based vs bio-based difference. For all analyzed sample, the majority of the OPAs  
369 identified were plasticizers. Indeed, during the formulation of plastic materials, OPAs are added  
370 at different proportions depending on the polymer and its usage. Most of the time, plasticizers  
371 are used extensively, reaching up 60% weight/weight (w/w) [59]. Here, the same number of  
372 additives detected in both PLA and PP materials (10 OPAs each) confirms that “bio-based” do  
373 not pledge “additives free”. Overall, lower amounts (Area (A.U)/mass ( $\mu\text{g}$ )) of additives were  
374 noticed in PLA samples (except BPB, DCHP and TBP), while, surprisingly, bisphenol B, a toxic  
375 endocrine disruptor analogue of the bisphenol A [10], has only been detected in PLA samples.  
376 Zimmerman et al., [19,60] tested the *in-vitro* toxicity and the chemical content of many bio-  
377 based materials (PLA items used as FCM). Their study detected priority concern compounds

378 (such as additives, oligomers, along with non-intentional added substances (NIAS)) and  
379 highlighted that bio-based materials are not safer than conventional plastics regarding their  
380 chemicals content.

381 Secondly, an inter-suppliers difference. This study underlined non-negligible differences on the  
382 identification of OPAs in two food-packaging products manufactured with the same polymer  
383 and that are intended to have the same use. That reflects the complexity to generalize and  
384 estimate the potentially dangerous additives in the studied material and thus the exposition of  
385 marine organisms and consumer, since the chemical composition of a material don't seem to  
386 be related to the polymer type but to each individual product. This match the conclusion of  
387 Zimmermann, et al. [61] that investigated the toxicity and the chemical composition of several  
388 plastic leachates made of conventional in PLA sample including FCM. They ended up with that  
389 same conclusion being that the contaminant content and the toxicity cannot be generalized  
390 based on the polymer type and are product specific. Plastic products formulation is different  
391 depending on the manufacturers, and may explain this difference in additives types and  
392 additives abundances between products made with the same polymer. It underlines the need  
393 of identifying the chemical content and the chemical safety of each material rather than  
394 polymers types, which is not a straightforward task but is inevitable considering the disparity  
395 in plastic formulation of the commercially available products. Moreover, the formulation's  
396 recipes remain unknown from the scientific community and hinder the evaluation of potential  
397 impacts and toxicity of plastic materials.

#### 398 **4. Conclusion**

399 This study was conducted with the aim to set up a suitable Py/TD-GC–HRMS method for  
400 simultaneous identification of various selected organic additives and to set up a high-resolution  
401 spectral database. The developed method showed its ability to detect the majority of the  
402 targeted additives. Nevertheless, this method may be limiting for certain OPAs harder to  
403 mobilize using GC. In complement, a liquid extraction analysis would provide more  
404 comprehensive results. Ten additives were detected in PP and PLA food contact materials.  
405 The results indicate the presence of multiples additives, some of them being toxic, in both  
406 petro-based and bio-based food packaging samples (e.g. plasticizers, phosphorous flames  
407 retardants). Moreover, significant differences were pointed out depending on either the  
408 polymer type or the suppliers. The present study highlights the complexity to assess and  
409 generalize the additive content and their amount in plastic packaging items due the variety and  
410 differences in plastic formulation.

## 411 **5. Declaration of competing interest**

412 Authors do not declare any competing interest.

## 413 **6. Funding**

414 Fleurine Akoueson is very grateful to the Hauts-de-France Regional Council and ANSES  
415 (French Agency for Food, Environmental and Occupational Health & Safety) for financially  
416 supporting her PhD scholarship. This work was also funded by the European Union (ERDF),  
417 the French government, the Hauts-de-France Regional Council and IFREMER, in the  
418 framework of the CPER MARCO and IDEAL projects. The funding of CPER MARCO project  
419 consists in the acquisition of the Py-GC-HRMS device and consumables used for the  
420 experiments.

## 421 **7. CRediT Roles**

422 **F. Akoueson**: Conceptualization, Investigation, Methodology, Formal analysis, Software,  
423 Validation, Visualization, Writing - original draft, Writing - review & editing. **C. Chbib**:  
424 Methodology Investigation, Writing - review & editing. **A. Brémard**: Investigation, Formal  
425 analysis. **S. Monchy**: Writing - review & editing. **I. Paul-Pont**: Writing - review & editing. **P.**  
426 **Doyen**: Writing - review & editing. **A. Dehaut**: Conceptualization, Funding acquisition,  
427 Supervision, Project administration, Writing - review & editing. **G. Duflos**: Resources,  
428 Conceptualization, Funding acquisition, Supervision, Project administration, Writing - review &  
429 editing.

## 430 **8. REFERENCES**

- 431 [1] K.J. Groh, T. Backhaus, B. Carney-Almroth, B. Geueke, P.A. Inostroza, A. Lennquist, H.A. Leslie,  
432 M. Maffini, D. Slunge, L. Trasande, A.M. Warhurst and J. Muncke, Overview of known plastic  
433 packaging-associated chemicals and their hazards, *Science of The Total Environment*, 651,  
434 (2019) 3253. <https://doi.org/10.1016/j.scitotenv.2018.10.015>
- 435 [2] PlasticEurope, Plastic Europe, the facts 2021, 2021, Plastic Europe, Available on:  
436 [https://plasticseurope.org/wp-content/uploads/2021/12/Plastics-the-Facts-2021-web-](https://plasticseurope.org/wp-content/uploads/2021/12/Plastics-the-Facts-2021-web-final.pdf)  
437 [final.pdf](https://plasticseurope.org/wp-content/uploads/2021/12/Plastics-the-Facts-2021-web-final.pdf)
- 438 [3] PlasticEurope, Plastic Europe, the facts 2020 2020, Plastic Europe, Available on:  
439 <https://plasticseurope.org/knowledge-hub/plastics-the-facts-2020/>
- 440 [4] O. Gerigny, M. Brun, C. Tomasino, M. Le Moigne, C. Lacroix, M. Kerambrun and F. Galgani,  
441 Évaluation du descripteur 10 « Déchets marins » en France métropolitaine. Rapport  
442 scientifique pour l'évaluation 2018 au titre de la DCSMM, 2018, Available on:  
443 [https://sextant.ifremer.fr/documentation/dcsmm/documents/Evaluation\\_2018/Rapport\\_Ev](https://sextant.ifremer.fr/documentation/dcsmm/documents/Evaluation_2018/Rapport_Evaluation_DCSMM_2018_D10>Ifremer_CEDRE.pdf)  
444 [aluation\\_DCSMM\\_2018\\_D10>Ifremer\\_CEDRE.pdf](https://sextant.ifremer.fr/documentation/dcsmm/documents/Evaluation_2018/Rapport_Evaluation_DCSMM_2018_D10>Ifremer_CEDRE.pdf)
- 445 [5] M. Kumar, H. Chen, S. Sarsaiya, S. Qin, H. Liu, M.K. Awasthi, S. Kumar, L. Singh, Z. Zhang, N.S.  
446 Bolan, A. Pandey, S. Varjani and M.J. Taherzadeh, Current research trends on micro- and nano-  
447 plastics as an emerging threat to global environment: A review, *J Hazard Mater*, 409, (2021)  
448 124967. <https://doi.org/10.1016/j.jhazmat.2020.124967>



- 449 [6] C. Campanale, C. Massarelli, I. Savino, V. Locaputo and V.F. Uricchio, A Detailed Review Study  
450 on Potential Effects of Microplastics and Additives of Concern on Human Health, *Int J Environ*  
451 *Res Public Health*, 17, (2020). <https://doi.org/10.3390/ijerph17041212>
- 452 [7] J.N. Hahladakis, C.A. Velis, R. Weber, E. Iacovidou and P. Purnell, An overview of chemical  
453 additives present in plastics: Migration, release, fate and environmental impact during their  
454 use, disposal and recycling, *J Hazard Mater*, 344, (2018)  
455 179. <https://doi.org/10.1016/j.jhazmat.2017.10.014>
- 456 [8] I.A. Sheikh and M.A. Beg, Structural characterization of potential endocrine disrupting activity  
457 of alternate plasticizers di-(2-ethylhexyl) adipate (DEHA), acetyl tributyl citrate (ATBC) and  
458 2,2,4-trimethyl 1,3-pentanediol diisobutyrate (TPIB) with human sex hormone-binding  
459 globulin, *Reprod Toxicol*, 83, (2019) 46. <https://doi.org/10.1016/j.reprotox.2018.11.003>
- 460 [9] D. Lithner, Å. Larsson and G. Dave, Environmental and health hazard ranking and assessment  
461 of plastic polymers based on chemical composition, *Science of The Total Environment*, 409,  
462 (2011) 3309. <https://doi.org/10.1016/j.scitotenv.2011.04.038>
- 463 [10] H. Serra, C. Beausoleil, R. Habert, C. Minier, N. Picard-Hagen and C. Michel, Evidence for  
464 Bisphenol B Endocrine Properties: Scientific and Regulatory Perspectives, *Environ Health*  
465 *Perspect*, 127, (2019) 106001. <https://doi.org/10.1289/EHP5200>
- 466 [11] S. Eladak, T. Grisin, D. Moison, M.J. Guerquin, T. N'Tumba-Byn, S. Pozzi-Gaudin, A. Benachi, G.  
467 Livera, V. Rouiller-Fabre and R. Habert, A new chapter in the bisphenol A story: bisphenol S and  
468 bisphenol F are not safe alternatives to this compound, *Fertil Steril*, 103, (2015)  
469 11. <https://doi.org/10.1016/j.fertnstert.2014.11.005>
- 470 [12] J. Muncke, A.M. Andersson, T. Backhaus, J.M. Boucher, B. Carney Almroth, A. Castillo Castillo,  
471 J. Chevrier, B.A. Demeneix, J.A. Emmanuel, J.B. Fini, D. Gee, B. Geueke, K. Groh, J.J. Heindel, J.  
472 Houlihan, C.D. Kassotis, C.F. Kwiatkowski, L.Y. Lefferts, M.V. Maffini, O.V. Martin, J.P. Myers,  
473 A. Nadal, C. Nerin, K.E. Pelch, S.R. Fernandez, R.M. Sargis, A.M. Soto, L. Trasande, L.N.  
474 Vandenberg, M. Wagner, C. Wu, R.T. Zoeller and M. Scheringer, Impacts of food contact  
475 chemicals on human health: a consensus statement, *Environ Health*, 19, (2020)  
476 25. <https://doi.org/10.1186/s12940-020-0572-5>
- 477 [13] A.A. Koelmans, E. Besseling and E.M. Foekema, Leaching of plastic additives to marine  
478 organisms, *Environ Pollut*, 187, (2014) 49. <http://dx.doi.org/10.1016/j.envpol.2013.12.013>
- 479 [14] C.M. Rochman, E. Hoh, T. Kurobe and S.J. Teh, Ingested plastic transfers hazardous chemicals  
480 to fish and induces hepatic stress, *Sci Rep*, 3, (2013) 3263. <https://doi.org/10.1038/srep03263>
- 481 [15] E. Fasano, F. Bono-Blay, T. Cirillo, P. Montuori and S. Lacorte, Migration of phthalates,  
482 alkylphenols, bisphenol A and di(2-ethylhexyl)adipate from food packaging, *Food Control*, 27,  
483 (2012) 132. <https://doi.org/10.1016/j.foodcont.2012.03.005>
- 484 [16] P. Schmid and F. Welle, Chemical Migration from Beverage Packaging Materials—A Review,  
485 *Beverages*, 6, (2020). <https://doi.org/10.3390/beverages6020037>
- 486 [17] J.A. Garde, R. Catala and R. Gavara, Global and specific migration of antioxidants from  
487 polypropylene films into food simulants, *J Food Prot*, 61, (1998)  
488 1000. <https://doi.org/10.4315/0362-028X-61.8.1000>
- 489 [18] Y. Tanaka, K. Takahashi, T. Sasao, T. Kirigaya, S. Hosoi, N. Nagaoka, T. Kawamura and H.  
490 Nakazawa, [Detection of tris(2,4-di-tert-butylphenyl)phosphite (Irgafos 168, antioxidant in  
491 plastics) and its oxide in commercial frozen vegetables], *Shokuhin Eiseigaku Zasshi*, 44, (2003)  
492 181. <https://doi.org/10.3358/shokueishi.44.181>
- 493 [19] L. Zimmermann, G. Dierkes, T.A. Ternes, C. Volker and M. Wagner, Benchmarking the in Vitro  
494 Toxicity and Chemical Composition of Plastic Consumer Products, *Environ Sci Technol*, 53,  
495 (2019) 11467. <https://doi.org/10.1021/acs.est.9b02293>
- 496 [20] L.S. Kato and C.A. Conte-Junior, Safety of Plastic Food Packaging: The Challenges about Non-  
497 Intentionally Added Substances (NIAS) Discovery, Identification and Risk Assessment, *Polymers*  
498 *(Basel)*, 13, (2021). <https://doi.org/10.3390/polym13132077>
- 499 [21] J.C.J. Bart, Polymer/additive analysis by flash pyrolysis techniques, *Journal of Analytical and*  
500 *Applied Pyrolysis*, 58-59, (2001) 3. [https://doi.org/10.1016/S0165-2370\(00\)00160-1](https://doi.org/10.1016/S0165-2370(00)00160-1)

- 501 [22] Y. Kudo, K. Obayashi, H. Yanagisawa, F. Maruyama, S. Fujimaki, H. Miyagawa and K.  
502 Nakagawa, Development of a screening method for phthalate esters in polymers using a  
503 quantitative database in combination with pyrolyzer/thermal desorption gas chromatography  
504 mass spectrometry, *Journal of Chromatography A*, 1602, (2019)  
505 441. <https://doi.org/10.1016/j.chroma.2019.06.014>
- 506 [23] F.C.-Y. Wang, Polymer additive analysis by pyrolysis–gas chromatography I. Plasticizers, *J.*  
507 *Chromatogr. A*, (2000) 12. [https://doi.org/10.1016/S0021-9673\(00\)00346-0](https://doi.org/10.1016/S0021-9673(00)00346-0)
- 508 [24] L. Hermabessiere and C.M. Rochman, Microwave-Assisted Extraction for Quantification of  
509 Microplastics Using Pyrolysis-Gas Chromatography/Mass Spectrometry, *Environ Toxicol Chem*,  
510 40, (2021) 2733. <https://doi.org/10.1002/etc.5179>
- 511 [25] F. Akoueson, C. Chbib, S. Monchy, I. Paul-Pont, P. Doyen, A. Dehaut and G. Duflos, Identification  
512 and quantification of plastic additives using pyrolysis-GC/MS: A review, *Sci Total Environ*, 773,  
513 (2021) 145073. <https://doi.org/10.1016/j.scitotenv.2021.145073>
- 514 [26] F. Maruyama, S. Fujimaki, Y. Sakamoto, Y. Kudo and H. Miyagawa, Screening of phthalates in  
515 polymer materials by pyrolysis GC/MS, *Anal Sci*, 31, (2015)  
516 3. <https://doi.org/10.2116/analsci.31.3>
- 517 [27] J.W. Kim, Y.M. Kim, H.M. Moon, A. Hosaka, C. Watanabe, N. Teramae, E.K. Choe and S.W.  
518 Myung, Comparative study of thermal desorption and solvent extraction-gas chromatography-  
519 mass spectrometric analysis for the quantification of phthalates in polymers, *J Chromatogr A*,  
520 1451, (2016) 33. <https://doi.org/10.1016/j.chroma.2016.05.014>
- 521 [28] M. Llana-Ruiz-Cabello, S. Pichardo, N.T. Jimenez-Morillo, F.J. Gonzalez-Vila, E. Guillamon, J.M.  
522 Bermudez, S. Aucejo, A.M. Camean and J.A. Gonzalez-Perez, Pyrolysis-gas chromatography-  
523 isotope ratio mass spectrometry for monitoring natural additives in polylactic acid active food  
524 packages, *J Chromatogr A*, 1525, (2017) 145. <https://doi.org/10.1016/j.chroma.2017.10.023>
- 525 [29] E.D. Okoffo, F. Ribeiro, J.W. O'Brien, S. O'Brien, B.J. Tschärke, M. Gallen, S. Samanipour, J.F.  
526 Mueller and K.V. Thomas, Identification and quantification of selected plastics in biosolids by  
527 pressurized liquid extraction combined with double-shot pyrolysis gas chromatography–mass  
528 spectrometry, *Science of The Total Environment*, 715, (2020)  
529 136924. <https://doi.org/10.1016/j.scitotenv.2020.136924>
- 530 [30] K.D. Jansson, C.P. Zawodny and T.P. Wampler, Determination of polymer additives using  
531 analytical pyrolysis, *Journal of Analytical and Applied Pyrolysis*, 79, (2007)  
532 353. <https://doi.org/10.1016/j.jaap.2006.12.009>
- 533 [31] M. Herrera, G. Matuschek and A. Kettrup, Fast identification of polymer additives by pyrolysis-  
534 gas chromatography/mass spectrometry, *Journal of Analytical and Applied Pyrolysis*, 70,  
535 (2003) 35. [https://doi.org/10.1016/S0165-2370\(02\)00078-5](https://doi.org/10.1016/S0165-2370(02)00078-5)
- 536 [32] International Electrotechnical Commission, IEC 62321-8 - Determination of certain substances  
537 in electrotechnical products - Part 8: Phthalates in polymers by gas chromatography-mass  
538 spectrometry (GC-MS), gas chromatography-mass spectrometry using a pyrolyzer/thermal  
539 desorption accessory (Py/TD-GC-MS), 2017, Available on:  
540 <http://mapnamagz.yabesh.ir/std/handle/yse/233172>
- 541 [33] Y.M. Kim, J.W. Kim, H.M. Moon, M.J. Lee, A. Hosaka, A. Watanabe, N. Teramae, Y.K. Park and  
542 S.W. Myung, Rapid Quantification of N-Methyl-2-pyrrolidone in Polymer Matrices by Thermal  
543 Desorption-GC/MS, *Anal Sci*, 33, (2017) 821. <https://doi.org/10.2116/analsci.33.821>
- 544 [34] H. Yanagisawa, Y. Kudo, K. Nakagawa, H. Miyagawa, F. Maruyama and S.  
545 Fujimaki, Simultaneous Screening of Major Flame Retardants and Plasticizers in Polymer  
546 Materials Using Pyrolyzer/Thermal Desorption Gas Chromatography Mass Spectrometry  
547 (Py/TD-GC-MS), *Molecules*, 23, (2018). <https://doi.org/10.3390/molecules23040728>
- 548 [35] P. Gimeno, S. Thomas, C. Bousquet, A.F. Maggio, C. Civade, C. Brenier and P.A.  
549 Bonnet, Identification and quantification of 14 phthalates and 5 non-phthalate plasticizers in  
550 PVC medical devices by GC-MS, *J Chromatogr B Analyt Technol Biomed Life Sci*, 949-950, (2014)  
551 99. <http://dx.doi.org/10.1016/j.jchromb.2013.12.037>

- 552 [36] J. Zhang, Phtalathes Analysis With Method GB 5009.271-2016 Using the Agilent 8890 GC and  
553 MSD with Agilent JetClean, 2019, Agilent Technologies Application Notes, Available on:  
554 [https://www.agilent.com/cs/library/applications/application-phthalates-8890-gc-jetclean-](https://www.agilent.com/cs/library/applications/application-phthalates-8890-gc-jetclean-5994-0656en-agilent.pdf)  
555 [5994-0656en-agilent.pdf](https://www.agilent.com/cs/library/applications/application-phthalates-8890-gc-jetclean-5994-0656en-agilent.pdf)
- 556 [37] C. Chattelier and F. Lestremau, Considérations sur certains aspects metrologiques liés à la  
557 mesure du 4-nonylphénol - Etat de l'art et perspectives – Rapport AQUAREF 2014, 2014,  
558 Available on: [https://professionnels.ofb.fr/sites/default/files/pdf/2014\\_109.pdf](https://professionnels.ofb.fr/sites/default/files/pdf/2014_109.pdf)
- 559 [38] M. Bolgar, J. Hubball, J. Groeger and S. Meronek, *Handbook for the chemical analysis of plastic*  
560 *and polymer additives*, Taylor & Francis Group, 2016, p.  
561 <https://doi.org/10.1201/9781420044881>
- 562 [39] R. Espejo, K. Valter, M. Simona, Y. Janin and P. Arrizabalaga, Determination of nineteen 4-  
563 alkylphenol endocrine disrupters in Geneva municipal sewage wastewater, *Journal of*  
564 *Chromatography A*, 976, (2002) 335. [https://doi.org/10.1016/S0021-9673\(02\)01229-3](https://doi.org/10.1016/S0021-9673(02)01229-3)
- 565 [40] M. Frederiksen, K. Vorkamp, R. Bossi, F. Rigét, M. Dam and B. Svensmark, Method development  
566 for simultaneous analysis of HBCD, TBBPA, and dimethyl-TBBPA in marine biota from  
567 Greenland and the Faroe Islands, *International Journal of Environmental Analytical Chemistry*,  
568 87, (2007) 1095. <https://doi.org/10.1080/03067310701451178>
- 569 [41] F. Akoueson, C. Chbib, A. Dehaut and G. Duflos, Organic Plastic Additives High-Resolution Mass  
570 Spectra., 2022, Zenodo, <https://doi.org/10.5281/zenodo.6599058>
- 571 [42] G. Castello, P. Moretti and S. Vezzani, Retention models for programmed gas chromatography,  
572 *J Chromatogr A*, 1216, (2009) 1607. <https://doi.org/10.1016/j.chroma.2008.11.049>
- 573 [43] P. Fields, S. Allen, Z. Korunic, A. McLaughlin and T. Stathers, *Standardised testing for*  
574 *diatomaceous earth*, at: 2003, 779.
- 575 [44] B. Kolahgar and E. Pfannkoch, The use of different PTV inlet liner types for trapping alkanes,  
576 aromatics and oxygenated compounds during thermal desorption; 2002; Available on:  
577 [https://www.gerstelus.com/wp-content/uploads/2017/07/technote-2002-03\\_cis-liner-](https://www.gerstelus.com/wp-content/uploads/2017/07/technote-2002-03_cis-liner-type.pdf)  
578 [type.pdf](https://www.gerstelus.com/wp-content/uploads/2017/07/technote-2002-03_cis-liner-type.pdf); Accessed on: 22/08/2022;
- 579 [45] B. Geueke, C.C. Wagner and J. Muncke, Food contact substances and chemicals of concern: a  
580 comparison of inventories, *Food Addit Contam Part A Chem Anal Control Expo Risk Assess*, 31,  
581 (2014) 1438. <https://doi.org/10.1080/19440049.2014.931600>
- 582 [46] R Core Team, R: The R Project for Statistical Computing; Available on: [https://www.r-](https://www.r-project.org/)  
583 [project.org/](https://www.r-project.org/); Accessed on: 2022-05-15;
- 584 [47] J. Fox, S. Weisberg, B. Price, D. Adler, D. Bates, G. Baud-Bovy, B. Bolker, S. Ellison, D. Firth, M.  
585 Friendly, A. Zeileis and R-Core, Companion to Applied Regression, *R package version*, (2022)
- 586 [48] F. De Mendiburu, Statistical Procedures for Agricultural Research, 2021, R package version,  
587 Available on: <https://cran.r-project.org/web/packages/agricolae/index.html>
- 588 [49] T. Pohlert, Calculate Pairwise Multiple Comparisons of Mean Rank Sums, 2021, R package  
589 version, Available on: <https://cran.r-project.org/web/packages/PMCMR/PMCMR.pdf>
- 590 [50] ECHA, Mapping exercise – Plastic additives initiative; 2019; Available on:  
591 <https://echa.europa.eu/fr/mapping-exercise-plastic-additives-initiative>; Accessed on:  
592 2020/07/02;
- 593 [51] J.H. Bridson, E.C. Gaugler, D.A. Smith, G.L. Northcott and S. Gaw, Leaching and extraction of  
594 additives from plastic pollution to inform environmental risk: A multidisciplinary review of  
595 analytical approaches, *J Hazard Mater*, 414, (2021)  
596 125571. <https://doi.org/10.1016/j.jhazmat.2021.125571>
- 597 [52] E. Fries, J.H. Dekiff, J. Willmeyer, M.-T. Nuelle, M. Ebert and D. Remy, Identification of polymer  
598 types and additives in marine microplastic particles using pyrolysis-GC/MS and scanning  
599 electron microscopy, *Environmental Science. Processes & Impacts*, 15, (2013)  
600 1949. <https://doi.org/10.1039/c3em00214d>
- 601 [53] F.C.-Y. Wang, Polymer additive analysis by pyrolysis–gas chromatography- II. Flame retardants,  
602 *Journal of Chromatography A*, 886, (2000) 225. [https://doi.org/10.1016/S0021-](https://doi.org/10.1016/S0021-9673(00)00538-0)  
603 [9673\(00\)00538-0](https://doi.org/10.1016/S0021-9673(00)00538-0)

- 604 [54] F.C.-Y. Wang, Polymer additive analysis by pyrolysis–gas chromatography - IV. Antioxydants,  
605 *Journal of Chromatography A*, 891, (2000) 325. [https://doi.org/10.1016/S0021-](https://doi.org/10.1016/S0021-9673(00)00647-6)  
606 [9673\(00\)00647-6](https://doi.org/10.1016/S0021-9673(00)00647-6)
- 607 [55] EFSA, Scientific Opinion of the Panel on food contact materials, enzymes, flavourings and  
608 processing aids (CEF) on 25th list of substances for food contact materials, in *EFSA Journal*,  
609 EFSA Journal edn., 2009, Chapter 1-20, p. 1196. <https://doi.org/10.2903/j.efsa.2009.1196>
- 610 [56] A. Venema and R.C. Boom-Van Geest, In-situ hydrolysis/methylation pyrolysis CGC for the  
611 characterization of polyaramides, *Journal of Microcolumn Separations*, 7, (1995) 337.  
612 <https://doi.org/10.1002/mcs.1220070406>
- 613 [57] I. IyigUndoGdu, G.A. UstUnda and Y. Duydu, Toxicological Evaluation of Bisphenol A and Its  
614 Analogues, *Turk J Pharm Sci*, 17, (2020) 457. <https://doi.org/10.4274/tjps.galenos.2019.58219>
- 615 [58] V. Becerra and J. Odermatt, Detection and quantification of traces of bisphenol A and bisphenol  
616 S in paper samples using analytical pyrolysis-GC/MS, *Analyst*, 137, (2012)  
617 2250. <https://doi.org/10.1039/C2AN15961A>
- 618 [59] E.L. Teuten, J.M. Saquing, D.R. Knappe, M.A. Barlaz, S. Jonsson, A. Bjorn, S.J. Rowland, R.C.  
619 Thompson, T.S. Galloway, R. Yamashita, D. Ochi, Y. Watanuki, C. Moore, P.H. Viet, T.S. Tana,  
620 M. Prudente, R. Boonyatumanond, M.P. Zakaria, K. Akkhavong, Y. Ogata, H. Hirai, S. Iwasa, K.  
621 Mizukawa, Y. Hagino, A. Imamura, M. Saha and H. Takada, Transport and release of chemicals  
622 from plastics to the environment and to wildlife, *Philos Trans R Soc Lond B Biol Sci*, 364, (2009)  
623 2027. <https://doi.org/10.1098/rstb.2008.0284>
- 624 [60] L. Zimmermann, A. Dombrowski, C. Volker and M. Wagner, Are bioplastics and plant-based  
625 materials safer than conventional plastics? In vitro toxicity and chemical composition, *Environ*  
626 *Int*, 145, (2020) 106066. <https://doi.org/10.1016/j.envint.2020.106066>
- 627 [61] L. Zimmermann, Z. Bartosova, K. Braun, J. Oehlmann, C. Volker and M. Wagner, Plastic Products  
628 Leach Chemicals That Induce In Vitro Toxicity under Realistic Use Conditions, *Environ Sci*  
629 *Technol*, 55, (2021) 11814. <https://doi.org/10.1021/acs.est.1c01103>




630

631 **Figure captions**

632 **Figure 1:** Chromatograms of the targeted plastic additives analyzed using GC-HRMS. (A)  
633 plasticizers additives; (B) flames retardants additives; and (C) antioxidants, UV stabilizers and  
634 plasticizers additives. The numbers refer to the name of the molecules mentioned in Table 1.

635 **Figure 2:** Peak areas of some selected OPAs depending on the thermal desorption  
636 temperature (250 °C, 350 °C and 450 °C, respectively in blue, yellow and grey). (A) plasticizers  
637 and flame-retardants; (B) antioxidants; and (C) antioxidants and UV-stabilizers. Values are  
638 expresses as mean± 95% confidence interval. Homogeneous groups are indicated by the  
639 same letter.(n=3).

640 **Figure 3:** Corrected peak areas of the OPAs found in the 4 food plastic packaging samples in  
641 polypropylene (PP) and polylactic acid (PLA). The letters A and B referred to the suppliers.  
642 Homogeneous groups are indicated by the same letter. (n=3).








Figure 3





643 Table 1: selected organic plastic additives (OPAs) with their corresponding suppliers and the quantity that allowed peaks detection using Py-GC-  
 644 HRMS analysis. The molecules tagged with “\*” represent the OPAs not detected (n.d) at the highest concentration available.

| Function          | N° | Additive                             | Abbreviation | CAS         | Suppliers    | Limit of detection (ng) (split 5) |              |              |              |
|-------------------|----|--------------------------------------|--------------|-------------|--------------|-----------------------------------|--------------|--------------|--------------|
| Plasticizers      | 1  | Dimethyl phthalate                   | DMP          | 131-11-3    | Restek       | 5                                 |              |              |              |
|                   | 2  | Diethyl phthalate                    | DEP          | 84-66-2     |              |                                   |              |              |              |
|                   | 3  | Di-allyl phthalate                   | DAIP         | 131-17-9    |              |                                   |              |              |              |
|                   | 4  | Diisobutyl phthalate                 | DIBP         | 84-69-5     |              |                                   |              |              |              |
|                   | 5  | Di-n-butyl phthalate                 | DBP          | 84-74-2     |              |                                   |              |              |              |
|                   | 6  | Tributyl Acetyl Citrate              | ATBC         | 77-90-7     | Accustandard |                                   |              |              |              |
|                   | 7  | Di-n-hexyl phthalate                 | DHP          | 84-75-3     | Restek       |                                   |              |              |              |
|                   | 8  | Benzyl butyl phthalate               | BBP          | 85-68-7     |              |                                   |              |              |              |
|                   | 9  | Bis-2-Ethylhexyl Adipate             | DEHA         | 103-23-1    | Accustandard |                                   |              |              |              |
|                   | 10 | Diisoheptyl phthalate                | DIHP         | 71888-89-6  | Accustandard |                                   | 25           |              |              |
|                   | 11 | Tri(2-ethylhexyl) phosphate          | TEHPA *      | 78-42-2     | Accustandard |                                   | n.d (at 100) |              |              |
|                   | 12 | Dicyclohexyl phthalate               | DCHP         | 84-61-7     | Restek       |                                   | 5            |              |              |
|                   | 13 | Bis(2-Ethylhexyl) phthalate          | DEHP         | 117-81-7    |              |                                   |              |              |              |
|                   | 14 | Diisononyl hexahydrophthalate        | DINCH        | 166412-78-8 | Accustandard |                                   |              |              |              |
|                   | 15 | Di-n-octyl phthalate                 | DIOP         | 117-84-0    | Restek       |                                   |              |              |              |
|                   | 16 | Diisononyl phthalate                 | DINP         | 68515-48-0  | Accustandard |                                   |              | 25           |              |
|                   | 17 | Di-nonyl phthalate                   | DNP          | 84-76-4     | Restek       |                                   |              | 5            |              |
|                   | 18 | Diisodecyl phthalate                 | DIDP         | 68515-49-1  | Accustandard |                                   |              | 25           |              |
| Flames retardants | 19 | Triethyl Phosphate                   | TEP          | 78-40-0     | LGC standard | 5                                 |              |              |              |
|                   | 20 | Tripropyl Phosphate                  | TPP          | 115-86-6    |              |                                   |              |              |              |
|                   | 21 | Tributyl Phosphate                   | TBP          | 126-73-8    |              |                                   |              |              |              |
|                   | 22 | 2,4,6-Tribromophenol                 | 2,4,6,TBP    | 118-79-6    |              |                                   |              | LGC standard |              |
|                   | 23 | Tris(2-Chloroethyl)Phosphate         | TCEP         | 115-96-8    |              |                                   |              | LGC standard |              |
|                   | 24 | Tris(2-Chloroisopropyl)Phosphate     | TCPP         | 13674-84-5  | LGC standard |                                   |              |              |              |
|                   | 25 | 2,4,4'-Tribromodiphenyl ether        | BDE-28       | 41318-75-6  | Accustandard |                                   |              | 10           |              |
|                   | 26 | Tris(1,3-Dichloro-2-Propyl)Phosphate | TDCPP        | 13674-87-8  | LGC standard |                                   |              | 5            |              |
|                   | 27 | Triphenyl Phosphate                  | TPhP         | 513-08-6    |              |                                   |              |              |              |
|                   | 28 | 2,2',4,4'-Tetrabromodiphenyl ether   | BDE-47       | 5436-43-1   | Accustandard |                                   |              |              | 10           |
|                   | 29 | Tricresyl Phosphate                  | TCP *        | 1330-78-5   | Accustandard |                                   |              |              | n.d (at 100) |
|                   | 30 | Tricresyl Phosphate - isomer         | TCrP *       | 78-30-8     |              |                                   | n.d (at 100) |              |              |
|                   | 31 | 2,2',4,4',6-Pentabromodiphenyl ether | BDE-100      | 60348-60-9  | Accustandard |                                   | 5            |              |              |

| Function   | N° | Additive   | Abbreviation          | CAS                | Suppliers    | Limit of detection (ng) (split 5) |
|--|----|--|-----------------------|--------------------|--------------|-----------------------------------|
|  | 32 | Tri-o-tolyl phosphate  | <b>TToP *</b>         | <b>78-30-8</b>     | Accustandard | n.d (at 100)                      |
|  | 33 | 2,2',4,4',5-Pentabromodiphenyl ether                           | <b>BDE-99</b>         | <b>189084-64-8</b> | Accustandard | 10                                |
|  | 34 | 2,2',4,4',5,5'-Hexabromodiphenyl ether                         | <b>BDE-153</b>        | <b>68631-49-2</b>  |              |                                   |
|  | 35 | 2,2',4,4',5,6'-Hexabromodiphenyl ether                         | <b>BDE-154</b>        | <b>207122-15-4</b> |              |                                   |
|  | 36 | 2,2',3,4,4',5',6-Heptabromodiphenyl ether                      | <b>BDE-183</b>        | <b>207122-16-5</b> |              |                                   |
|  | 37 | 1,2-Bis (2,4,6 Tribromophenoxy) ethane                         | <b>BTBPE</b>          | <b>37853-59-1</b>  |              |                                   |
| <b>Antioxidants</b>                              | 38 | 6,6'-di-tert-butyl-2,2'-thiodi-p-cresol                        | <b>Irganox® 1081</b>  | <b>90-66-4</b>     | Accustandard | 10                                |
|  | 39 | Butylated hydroxytoluene                                       | <b>BHT</b>            | <b>128-37-0</b>    | Accustandard | 5                                 |
|  | 40 | 3,5-di-tert-butyl-4-hydroxyhydrocinnamic acid, octadecyl ester | <b>Irganox® 1076</b>  | <b>2082-79-3</b>   | Accustandard | 100                               |
|  | 41 | 6,6'-ditert-butyl-4,4'-thiodin-m-cresol                        | <b>Lowinox® TBM-6</b> | <b>96-69-5</b>     |              | 50                                |
| <b>UV stabilizers</b>                            | 42 | 2,2-dihydroxy-4,4-dimethoxybenzophenone                        | <b>Uvinul® 3049 *</b> | <b>131-54-4</b>    | Accustandard | n.d. (at 100)                     |
|  | 43 | 2-t-Butyl-6(5-chloro-2H-benzotriazol-2-yl)-4-methylphenol      | <b>UV-326</b>         | <b>3896-11-5-</b>  |              | 5                                 |
|  | 44 | 2-(2H-Benzotriazol-2-yl)-4,6-di-tert-pentylphenol              | <b>UV-328</b>         | <b>25973-55-1</b>  |              |                                   |
|  | 45 | 2,4-Di-tert-butyl-6-(5-chloro-2H-benzotriazol-2-yl)phenol      | <b>UV-327</b>         | <b>3864-99-1</b>   |              |                                   |
|  | 46 | 2-hydroxy-4-octyloxybenzophenone                               | <b>Uvinul 3008 *</b>  | <b>1843-05-6</b>   |              | n.d. (at 100)                     |
| <b>Antioxidants - plasticizers - stabilizers</b> | 47 | 4-Tert-Octylphenol   | <b>4-t-OP</b>         | <b>140-66-9</b>    | Accustandard | 5                                 |
|  | 48 | Nonylphenol  | <b>NPs</b>            | <b>84852-15-3</b>  | Accustandard | 10                                |
|  | 49 | 4-nonylphenol  | <b>4-NP</b>           | <b>104-40-5</b>    | Accustandard |                                   |
|  | 50 | Nonylphenol Monoethoxylate                                     | <b>NP1EO</b>          | <b>27986-36-3</b>  | Accustandard |                                   |
|  | 51 | Bisphenol F  | <b>BPF</b>            | <b>620-92-8</b>    | Accustandard | 5                                 |
|  | 52 | 4-Nonylphenol Monoethoxylate                                   | <b>4-NP1EO</b>        | <b>104-35-8</b>    | Accustandard | 10                                |
|  | 53 | Bisphenol A  | <b>BPA</b>            | <b>80-05-7</b>     | Accustandard | 5                                 |
|  | 54 | Bisphenol B  | <b>BPB</b>            | <b>77-40-7</b>     |              |                                   |
|  | 55 | Nonylphenol diethoxylate                                       | <b>NP2EO *</b>        | <b>N/A</b>         | Accustandard | n.d. (at 10)                      |
|  | 56 | Bisphenol S  | <b>BPS *</b>          | <b>80-09-1</b>     | Accustandard | n.d (at 100)                      |

645

646 Table 2: GC-HRMS analytical conditions

| Parameters                 | Settings                                  |
|----------------------------|---|
| <b>GC</b>                  |   |
| Injection port temperature | 300°C                                     |
| Column temperature (°C)    | 80°C (0.5 min) → 10°C/min → 330°C (3 min) |
| Injection mode             | Split (split ratio 1:5)                   |
| Carrier gas                | Helium                                    |
| Carrier Gas flow           | 1 mL.min <sup>-1</sup>                    |
| Transfer line temperature  | 300°C                                     |
| <b>MS</b>                  |   |
| Ion source temperature     | 300°C                                     |
| Ionization method          | Electron Ionization (EI), 70 eV.          |
| Analytical mode            | Full Scan mode (FS)                       |
| Scan Range                 | 30.00000 – 750.00000 m/z                  |
| Solvent Delay              | 5 min                                     |

647

648 **Table 3:** Kováts retention indices (RI) values and quantitative and qualitative ions for the  
 649 selected plastic additives molecules. n.d. = not detected.

| Function                                      | Molecules      | RI          | Quantitative ions | Qualitative ions      |
|---|----------------|-------------|-------------------|-----------------------|
| Plasticizers                                  | DMP            | 1469        | 163.03888         | 133.02846             |
|   | DEP            | 1608        | 149.02333         | 177.05455             |
|   | DAIP           | 1756        | 149.0231          | 189.05432             |
|   | DIBP           | 1872        | 149.02333         | 161.05917             |
|   | DBP            | 1972        | 149.02333         | 160.08830             |
|   | ATBC           | 2160        | 129.01845         | 259.15451 ; 137.99502 |
|   | DHP            | 2374        | 149.02333         | 234.12006             |
|   | BBP            | 2401        | 149.02333         | 206.09376             |
|   | DEHA           | 2404        | 111.04393         | 101.05960 ; 129.05444 |
|   | DIHP           | 2380 - 2509 | 149.02333         | 265.14313 ; 99.11672  |
|   | TEHPA          | 2477        | 98.98428          | 113.13261             |
|   | DCHP           | 2558        | 149.02333         | 249.11210             |
|   | DEHP           | 2570        | 149.02333         | 279.15909             |
|   | DINCH          | 2630 – 2870 | 155.07002         | 109.06474 ; 127.07529 |
|   | DIOP           | 2756        | 149.02333         | 261.14876             |
|   | DINP           | 2731 – 2867 | 149.02323         | 71.08553              |
|   | DNP            | 2936        | 149.02333         | 293.17474             |
| DIDP  | 2876 - 3010    | 149.02332   | n.d.              |                       |
| Flames retardants                             | TEP            | 1124        | 98.98432          | 127.01567 ; 109.00512 |
|   | TPP            | 1383        | 98.98432          | 141.03133             |
|   | TBP            | 1650        | 98.98432          | 211.10963             |
|   | 2,4,6,TBP      | 1656        | 331.76931         | 222.85798             |
|   | TCEP           | 1769        | 142.96611         | 222.96906 ; 248.98479 |
|   | TCPP           | 1808        | 125.00023         | 201.00710 ; 156.98178 |
|   | BDE-28         | 2283        | 245.96733         | 139.05414             |
|   | TDCPP          | 2354        | 98.98432          | 190.94284 ; 154.96611 |
|   | TPhP           | 2421        | 325.06296         | 233.03670 ; 215.02596 |
|   | BDE-47         | 2550        | 325.87582         | 323.87792             |
|   | TCP            | 2676        | 165.07015         | 243.05746 ; 368.11771 |
|   | TCrP           | 2707        |                   |                       |
|   | BDE-100        | 2763        | 403.78613         | 563.62171             |
|   | TTOP           | 2738        | 165.07015         | 243.05746 ; 368.11771 |
|   | BDE-99         | 2823        | 405.78452         | 563.62171             |
|   | BDE-153        | 3000        | 483.69452         | 242.84599 ; 643.53162 |
|   | BDE-154        | 3092        |                   |                       |
| BDE-183                                       | 3375           | 563.60352   | 281.80109         |                       |
| BTBPE   | 3453           | 356.79514   | 277.87637         |                       |
| Antioxidants                                  | Irganox® 1081  | 1366        | 149.09674         | 164.12027             |
|   | BHT            | 1520        | 205.158752        | 206.16206             |
|   | Irganox® 1076  | 3611        | 147.08102         | 263.16534 ; 217.15947 |
|   | Lowinox® TBM-6 | 2795        | 358.19647         | 259.07922             |
| UV stabilizers                                | Uvinul® 3049   | 2503        | 151.03903         | 257.08099 ; 124.05204 |
|   | UV-326         | 2575        | 300.08972         | 272.05853             |
|   | UV-328         | 2681        | 322.19113         | 351.23003             |
|   | UV-327         | 2692        | 342.13678         | 357.16034 ; 286.07428 |
|   | Uvinul 3008    | 2790        | 213.05461         | 325.17978 ; 137.02342 |
| Antioxidants;<br>Plasticizers;<br>Stabilizers | 4-t-OP         | 1611        | 135.08041         | 206.16644             |
|   | NPs            | 1700 – 1770 | 135.08041         | 197.13233 ; 212.15582 |
|   | 4-NP           | 1873        | 107.04916         | 221.18555             |
|   | NP1EO          | 1969 - 2030 | 179.10652         | 193.12220             |
|   | BPF            | 2104        | 200.08276         | n.d.                  |
|   | 4-NP1EO        | 2146        | 107.04916         | 151.07530 ; 264.20834 |
|   | BPA            | 2221        | 213.09113         | 228.11452             |
| BPB   | 2321           | 213.09065   | 227.10631         |                       |

| <b>Function</b> | <b>Molecules</b> | <b>RI</b>   | <b>Quantitative ions</b> | <b>Qualitative ions</b> |
|-----------------|------------------|-------------|--------------------------|-------------------------|
|                 | <b>NP2EO</b>     | 2027 - 2321 | 107.04916                | 151.07530 ; 264.20834   |
|                 | <b>BPS</b>       | 2765.2777   | 141.00046                | 110.03627               |

650

651

652 Table 4: average extraction efficiency percentages and, in parenthesis, standard deviations of  
 653 the method following the 1st pyrolysis step, after analyzing the additive content of PLA and PP  
 654 samples (n=3). Each sample was successively analyzed 3 times at 350 °C. Bold printed  
 655 characters represent average extraction efficiencies >80%. n.d = not detected.

|             | PLA-A          | PLA-B          | PP-A           | PP-B           |
|-------------|----------------|----------------|----------------|----------------|
| <b>DMP</b>  | n.d            | 33 (56)        | 79 (10)        | 48 (9)         |
| <b>DEP</b>  | <b>90 (17)</b> | <b>95 (9)</b>  | 70 (4)         | 65 (9)         |
| <b>DBP</b>  | <b>94 (10)</b> | <b>99 (3)</b>  | <b>100 (0)</b> | <b>100 (0)</b> |
| <b>DIBP</b> | <b>98 (3)</b>  | <b>98 (3)</b>  | <b>100 (0)</b> | <b>100 (0)</b> |
| <b>DCHP</b> | <b>100 (0)</b> | <b>100 (0)</b> | <b>100 (0)</b> | n.d            |
| <b>DINP</b> | n.d            | n.d            | <b>87 (2)</b>  | <b>88 (4)</b>  |
| <b>DEHA</b> | n.d            | <b>90 (18)</b> | 59 (36)        | <b>93 (2)</b>  |
| <b>TBP</b>  | 59 (51)        | <b>93 (6)</b>  | 23 (40)        | n.d            |
| <b>TCPP</b> | 62 (54)        | 76 (9)         | 65 (5)         | 79 (11)        |
| <b>BPB</b>  | <b>100 (0)</b> | <b>93 (2)</b>  | n.d            | n.d            |

656



[Click here to access/download](#)

**Supplementary Material**

**Supplementary datas - Akoueson et al.pdf**

