

RAPID, ON SITE, DIAGNOSIS OF FMD AND SAFE AND COST-EFFECTIVE SHIPMENT OF SAMPLES USING LATERAL FLOW DEVICES FOR LABORATORY DIAGNOSTIC

Labib Bakkali Kassimi, Graham J. Belsham, Abdulnaci Bulut, Kamila Gorna, Claude Hamers, Pascal Hudelet, Syed Jamal, Eve Laloy, Anthony Relmy, Aurore Romey, et al.

▶ To cite this version:

Labib Bakkali Kassimi, Graham J. Belsham, Abdulnaci Bulut, Kamila Gorna, Claude Hamers, et al.. RAPID, ON SITE, DIAGNOSIS OF FMD AND SAFE AND COST-EFFECTIVE SHIPMENT OF SAMPLES USING LATERAL FLOW DEVICES FOR LABORATORY DIAGNOSTIC. 5th Congress of the European Association of Veterinary Laboratory Diagnosticians, Oct 2018, Bruxelles, France. anses-04441008

HAL Id: anses-04441008 https://anses.hal.science/anses-04441008

Submitted on 6 Feb 2024

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



RAPID, ON SITE, DIAGNOSIS OF FMD AND SAFE AND COST-EFFECTIVE SHIPMENT OF SAMPLES USING LATERAL FLOW DEVICES FOR LABORATORY DIAGNOSTIC



Bakkali Kassimi L¹, Belsham, G.J.², Bulut A.N.³, Gorna, K.¹, Hamers, C.⁴, Hudelet, P.⁴, Jamal, S.⁵, Laloy, E.¹, Relmy, A.¹, Romey, A.¹, Ularamu H.G.⁶, Zientara, S.¹ and Blaise-Boisseau S.¹

¹ Laboratoire de Santé Animale de Maisons- Alfort, Laboratoire de référence Nationale et OIE pour la Fièvre Aphteuse, UMR Virologie 1161, Université Paris-Est, Anses, Maisons-Alfort, France; ² DTU National Veterinary Institute, Technical University of Denmark, Lindholm, Denmark (DTU-Vet); ³ SAP/FMD Institute, Dumlupinar Bulvard,35, 06510, Ankara, Turkey; ⁴ The Veterinary Public Health Center, Boehringer Ingelheim Animal Health, 29 Avenue Tony Garnier, 69007 Lyon France; ⁵ Department of Biotechnology, University of Malakand (UM), Chakdara, Pakistan; ⁶ FMD Research Centre, Nat. Vet. Res. Inst. (NVRI), PMB 01 Vom. Nigeria

Email: aurore.romey@anses.fr

Foot-and-mouth disease: FMD

- Highly contagious transboundary disease of domestic and wild cloven-hoofed animals
- One of the most economically devastating diseases of livestock
- Characterized by vesicles, with subsequent erosions in the mouth, also on the muzzle, feet, or teats
- High morbidity, low mortality in adults










Foot-and-mouth disease virus: FMDV


- Causative agent = FMDV
- *Picornaviridae*, genus *Aphtovirus*, genome: positive-sense, single-stranded RNA
- 7 serotypes (O/A/C/Asia-1/SAT1,2,3) sub-divided into multiple topotypes and genotypes
- 7 pools defined according to the world distribution of FMDV






Identification of circulating strains is important for efficient FMD control. However, due to shipping requirements of potentially infected samples, the cost of sample submission to reference laboratories remains a major obstacle.






A cost-effective and safe method for shipment of samples from FMD-suspected cases, based on the inactivation of FMDV on lateral flow devices (LFDs) has been developed and validated in the laboratory which allows subsequent detection and typing of FMDV by RT-PCR and virus rescue using RNA transfection (*Romey et al., TBED, 2017*).


Evaluation of this protocol on freshly collected clinical samples through collaborations in endemic countries.






Optimization of RNA chemical transfection to rescue infectious FMDV from positive inactivated LFDs

- Sample loading on LFD
- Inactivation of FMDV on LFD
- RNA extraction from LFD
- RNA transfection:
- 5 transfection reagents
- 8 conditions per reagent
- 3 cells lines



Mock



CPE on ZZ-R-127 48h post-transfection



- 2 transfection reagents
- 1 condition per reagent
- 1 cell line
- ... in progress

This project will allow to ensure that the use of LFD is a safe way for room-temperature dry-transport of FMDV samples from endemic areas. Once validated, this method will substantially decrease the cost of shipment and help in increasing sample submission from affected countries.











